

Seating Series

Primacare Motion Recliner™



**GLOBALcare**  
Quality of Life Products™

*Adjustable comfort with a mobility option.*



Upright position.



Reclined position (partial).  
Allows user to adjust the degree of the recline to best suit the comfort level.


**GC3609** Primacare Motion Recliner™ with 1/3 front solid wood armcap and 2/3 upholstered back armcap shown with wood leg option and optional IV pole (IL).

Primacare™ Motion Recliners are easy to operate and suit a wide range of healthcare applications. The upholstery details and shapes minimize dirt traps and facilitate easy cleaning and disinfecting in any position. The tailored appearance is compatible with a wide range of interior and architectural styles.

Primacare™ Motion Recliners incorporate high density “UltraCell” bio foam in the seat, back and 3 piece leg supports for superior comfort and longevity. Weight capacity is warranted up to 350 lbs.

Above: GC3609 shown in Mayer, Finish Line Greige (WC761-007), seat in Maharam, Ledger Loden (463770-017), armcaps and legs in Avant Maple (AVM).

*Features and Benefits*



Stationary Base Option - Decorative tapered wood legs (WL). Available in a variety of GLOBALcare wood finishes.



Mobile Base Option - Institutional quality locking casters (MRC). Not intended for patient transfer purposes.





Trendelenburg Handle - Can be easily operated by a single caregiver.

**GC3608TR** shown in optional Trendelenburg Position. Allows the caregiver to recline the chair to position the patient's feet, higher than the head. Shown in Mayer, Finish Line Greige (WC761-007), seat in Maharam, Ledger Loden (463770-017), with optional armcaps in Toast (TOA).


*Features and Benefits*



**GC3607** Optional contoured headrest provides exceptional comfort.



Chrome push bar is standard with the mobile base option (MRC).



Removable upholstery covers are easy to clean and replace.



Removable, zippered seat cushions allow for easy cleaning and cushion replacement.




## Armcaps




### Solid Wood Full Armcap

Available in a variety of GLOBALcare wood finishes. Field replaceable.



### Self-Skinned Urethane Armcap

Available in Black (BLK), Platinum (PLT) or Toast (TOA). Field replaceable.



### Solid Wood Half Armcap

Available in a variety of GLOBALcare wood finishes. Field replaceable.

#### DESCRIPTION

DIMENSIONS	DIMENSIONS			SEAT
	W	D	H	HEIGHT
	IN	IN	IN	IN
	mm	mm	mm	mm



<b>GC3608</b> Primacare motion recliner with solid wood full span armcaps or self skinned urethane (SSU) armcaps.	30	39	46.5	19
<b>GC3608MRC</b> Motion recliner with casters.	762	991	1181	483
<b>GC3608WL</b> Motion recliner with wood legs.				



<b>GC3607</b> Primacare motion recliner headrest.	14	4	9.5	N/A
	356	102	241	N/A

#### DESCRIPTION

DIMENSIONS	DIMENSIONS			SEAT
	W	D	H	HEIGHT
	IN	IN	IN	IN
	mm	mm	mm	mm



<b>GC3609</b> Primacare motion recliner with 1/3 front solid wood armcaps.	30	39	46.5	19
<b>GC3609MRC</b> Motion recliner with casters.	762	991	1181	483
<b>GC3609WL</b> Motion recliner with wood legs.				

Cover: GC3608 shown in Mayer, Finish Line Greige (WC761-007), seat in Maharam, Ledger Loden (463770-017), with optional SSU armcaps in Toast (TOA).

This brochure provides an overview of the key features and benefits of the GLOBALcare product. For a full list of products available, including specifications, fabrics and finishes, please visit [www.thinkglobalcare.com](http://www.thinkglobalcare.com). Please refer to our GLOBALcare binder for pricing inquiries or contact your local GLOBALcare sales representative.

**Global - Canada**  
The Global Group, 1350 Flint Rd., Downsview ON Canada M3J 2J7  
Sales & Marketing: Tel (1-877) 446-2251  
Customer Service: Fax (800) 361-3182  
Government Customer Service: Fax (416) 739-6319

**Global - U.S.A.**  
17 West Stow Road P.O. Box 562  
Marlton New Jersey U.S.A. 08053  
Tel (856) 596-3390 (800) 220-1900  
Fax (856) 596-5684

**Global - International**  
560 Supertest Road Downsview  
Ontario Canada M3J 2M6  
Tel (416) 661-3660 (800) 668-5870  
Fax (416) 667-1821

**GGI Office Furniture (UK) Ltd.**  
Global Way Darwen  
Lancashire UK BB3 0RW  
Tel +44(0) 1254 778500  
Fax +44(0) 1254 778519

Visit us on the World Wide Web at [www.thinkglobalcare.com](http://www.thinkglobalcare.com)

Offices throughout Canada, United States, United Kingdom, France, Brazil, Australia and Italy. Distribution around the world. A member of the Global Group.

