



MOTION RECLINER
primacare

healthcare

global

Recline effortlessly, transfer safely



Primacare™ Motion Recliners are easy to operate and suited to a wide range of healthcare applications from acute to ambulatory care and long-term care. The series includes models with manual locking casters or four wheel central locking casters accessible from both sides of the chair. Select stationary and mobile models are available with folding arms to safely support patient transfer.

Relax and recline

3 position recliner with optional reverse recline



First position



Second position



Third position



Reverse recline position (optional)

Primacare Motion Recliners follow the patient's intuitive response to reach for the arms when they want to come to a seated position. Additional ergonomic considerations include grip assist arms, waterfall front edge seat and built-in lumbar support. Primacare Motion Recliners are tested and rated for an active weight capacity of 400 lbs. Backed by Global's Limited Lifetime Warranty.

Additional information can be found on the Global [website](#).

Ease and accessibility

Tablet movements



Parallel to chair arm



45° to chair arm



90° to chair arm

An optional Black phenolic folding tablet with a spill collector groove can be specified on the right, left or both sides of a fixed arm. When locked in the upright position, the tablet rotates 360° with a controlled resistance, without obstructing the patient when sitting or exiting the chair.

Fold-down arm



Open fold-down arm



Closed fold-down arm



Caregiver activated fold-down arm release

Fold-down arms provide ease of access to assist with patient transfer. The auto-retracting fabric panel enhances patient protection and provides an easily cleanable surface. The arm automatically locks when placed in the upright position.

Fixed or folding arms define each model



Motion recliner with urethane armcaps on wood legs



Motion recliner with left fold-down arm, wood armcaps, right tablet, reverse recline and foot tray, on manual locking casters



Motion recliner with right fold-down arm, solid surface armcaps, left tablet, reverse recline and foot tray, on central locking casters



Motion recliner with right and left fold-down arms, solid surface armcaps, reverse recline and foot tray, on central locking casters

Models are available with fixed wood legs, fixed locking front and swivel rear casters or four swivel, locking casters. Select models feature fold-down arms (right, left or both) and an optional extendable foot tray to support the user's feet when seated.

Seating shown in Maharam Elm, Fencepost with LDI California, Misty Gray seat; Shadow urethane, WC Cream wood, Beige Tempest solid surface armcaps and Night wood legs.

A closer look



The Primacare Motion Recliner can be specified to best suit the needs of the patient, caregiver, and the intended application. Sealed seams are available on approved coated textiles for the back, arms and seat. Two textile specifications are available, as shown above.

It's in the details

Back



Removable cloth covering



Adjustable contoured headrest, ordered separately



Chrome push bar (optional)

Caregiver activated



Reverse recline foot lever (optional)

Storage



Swivel accessory hook (optional)

Features



Clean out on sides of seat cushions



Removable upholstery covers



Seat and back deck are constructed using a heavy-gauge, no-sag spring system supported by a painted steel frame


Armcaps



Urethane



Wood (optional)



Solid surface (optional)

Feet and casters



Wood



Manual locking casters (optional)



Central locking casters with foot paddle

Please visit globalfurnituregroup.com for additional product information including industry and environmental certifications.

Cover and above: Seating shown in Maharam Elm, Fencepost back and arms with LDI California, Misty Grey seat, Beige Tempest solid surface armcaps.



U.S.A.
Global Furniture Group
17 West Stow Road P.O. Box 562
Marlton New Jersey U.S.A. 08053
Tel (856) 596-2390 (800) 220-1900

Canada
Global Furniture Group
1350 Flint Rd., Toronto Ontario Canada M3J 2J7
Sales & Marketing: Tel (1-877) 446-2251

International
Global Furniture Group
560 Supertest Road, Toronto
Ontario Canada M3J 2M6
Tel (416) 661-3660 (800) 668-5870

Offices throughout Canada and the United States. Distribution around the world. A member of Global Furniture Group. © 2025 ®™ Trademarks of Global Furniture Group. These products may be covered by one or more issued patents or pending patent applications in Canada, the U.S. and elsewhere. All intellectual property rights in applicable trademarks, patents and industrial designs are strictly reserved. Made in Canada. Specifications are correct at time of printing. Printed in Canada.
Creative / Global Design Center 271025 24.0049